



Test Bulletin

February 2019

ACL Laboratories Offers Comprehensive Drug Screen – Pain Management Profile 1 (Test Order Code QPNMG1) Effective Wednesday, February 20, 2019, ACL Laboratories will offer Pain Management Profile 1 (Test Order Code QPNMG1), a comprehensive drug screen. With this new offering, the following pain management tests will be deactivated:

- Drug Abuse Panel – Drug Management (Test Order Code DMABUS)
- Pain Management Drug Panel (Test Order Code DRGPM2)
- Pain Management Drug Screen with Interpretation (Test Order Code DRUGPM).

Pain Management Profile 1 (Test Order Code QPNMG1), a comprehensive drug screen, is intended for use in pain management patients and those requiring routine drug monitoring. The panel utilizes immunoassay screening methods and definitive LC/MS methods to test for common prescribed and abused drugs. All positive immunoassay results auto-reflex to definitive LC/MS testing, and all positive results are quantitatively reported. Patient medication information must be provided to receive an interpretative report. Results will indicate whether findings are consistent or inconsistent with supplied medication information.

The specimen requirements are as follows:

Volume	30.0 mL Random Urine
Collection	Must be collected in a sterile urine container. Do not transfer into an aliquot tube
Transport	Refrigerated

The following drugs/classes of drugs will be screened and ALL positive results will auto-reflex to definitive LC/MS confirmation:

	Positive Cut-off
Amphetamine/Metamphetamine	500 ng/mL
Barbiturates	300 ng/mL
Buprenorphine/Naloxone	5 ng/mL
Cocaine Metabolite (Benzoylecgonine)	150 ng/mL
Fentanyl	0.5 ng/mL
Heroin Metabolite (6-Acetylmorphine)	10 ng/mL
Marijuana Metabolite (THC)	20 ng/mL
MDMA/MDA	500 ng/mL
Methadone	100 ng/mL
Phencyclidine	25 ng/mL
Propoxyphene	300 ng/mL

The following drugs will be tested definitively by LC/MS and all values above the cut-off concentration quantitatively reported:

	Positive Cut-off
Gabapentin	1,000 ng/mL
Meperidine & Normeperidine	100 ng/mL
Meprobamate (Carisoprodol metabolite)	1,000 ng/mL
Pregabalin	1,000 ng/mL
Tapentadol & Nortapentadol	50 ng/mL
Tramadol & Desmethyltramadol	100 ng/mL
Benzodiazepines, Quantitation	
Alphahydroxyalprazolam	25 ng/mL
Alphahydroxymidazolam	50 ng/mL
Alphahydroxytriazolam	50 ng/mL
Aminoclonazepam	25 ng/mL
Hydroxyethylflurazepam	50 ng/mL
Lorazepam	50 ng/mL
Nordiazepam	50 ng/mL
Oxazepam	50 ng/mL
Temazepam	50 ng/mL
Opiates, Expanded Quantitation	
Codeine	50 ng/mL
Hydrocodone	50 ng/mL
Hydromorphone	50 ng/mL
Morphine	50 ng/mL
Norhydrocodone	50 ng/mL
Noroxycodone	50 ng/mL
Oxycodone	50 ng/mL
Oxymorphone	50 ng/mL

For additional information regarding the Pain Management Profile 1 (Test Order Code **QPNMG1**) test, please visit ACL Laboratories Directory of Services at <http://www.acllaboratories.com/test-catalog/> or contact ACL Client Services at 1.800.877.7016.