









Microbiology Specimen Collection

| | | |
|---|---|---|
| ESwab – White Cap <i>For anaerobic, aerobic, MRSA/MSSA PCR, H1N1 and viruses by PCR – Ambient, SwabOne™ Panels</i> |  | |
| ESwab – Green Cap or Blue Cap <i>For nasopharyngeal or urethral sites for culture; Bordetella pertussis – Ambient</i> |  | |
| Ova & Parasite Transport SAF Media <i>For O&P, Cyclospora, Cystoisospora, Microsporidia – Ambient</i> |  | |
| Stool Culture Cary – Blair Media <i>For Salmonella, Shigella, Aeromonas, Plesiomonas, Campylobacter, Shiga Toxin, Vibrio, Yersinia – Ambient</i> <i>Giardia, Cryptosporidium, E. histolytica by PCR – Ambient</i> |  | |
| Aptima Unisex Swab – Female <i>For Chlamydia/GC and/or Trichomonas vaginalis Nucleic Acid Amplification – Ambient or Refrigerated</i> | Aptima Unisex Swab – Male <i>For Chlamydia/GC only – Ambient or Refrigerated</i> |  |
| Aptima Urine – Female or Male <i>For Chlamydia/GC and/or Trichomonas vaginalis Nucleic Acid Amplification – Ambient or Refrigerated</i> | |  |
| Aptima Vaginal Swab <i>For Chlamydia/GC Nucleic Acid Amplification, Throat or Rectal (Female or Male)</i> <i>For Chlamydia/GC Trichomonas vaginalis, Vaginal (Female only) – Ambient or Refrigerated</i> | |  |
| Nasopharyngeal Wire Swab <i>For Rapid RSV and Rapid Influenza only – Ambient</i> | |  |
| UTM Transport Media <i>Media for Mycoplasma, Ureaplasma, C. trachomatis Culture, and Mumps PCR or MERS CoV – Refrigerated</i> Rapid RSV – Ambient Rapid Influenza – Ambient, SwabOne™ Panels, Refrigerated | |  |
| Anaerobic and Aerobic Blood Culture Bottles Ambient | |  |
| Pediatric Blood Culture Bottle Ambient | |  |
| Pinworm Paddle Ambient | |  |

Urine Specimen Collection

| | |
|---|---|
| 4 mL Urine Collection Tube with Straw Submit 4 mL gray top when requesting URINE CULTURE. Ambient |  |
| 8 mL Urine Collection Tube with Straw Submit 8 mL red/yellow top when requesting URINALYSIS. – Ambient |  |
| Urine Transport Tube Refrigerated |  |
| 4 mL AND 8 mL Urine Collection Tubes with Straw Submit both 4 mL gray and 8 mL red/yellow top tubes when requesting urinalysis AND Urine Culture or Urinalysis with Culture, if indicated. – Ambient |  |

Clinical/Pathology Specimen Collection

| | |
|---|---|
| Sterile Container <i>For C. difficile Toxin Assay, Stool for WBCs, Occult Blood – Card Test, Rotavirus, Cytology Body Fluids</i> Refrigerated |  |
| FOBT – Chek <i>For Occult Blood – Immunochemical</i> Ambient |  |
| Quantiferon TB <i>Special handling required. Shake tubes firmly 10 times. Refer to Directory of Services or contact Client Services.</i> – Ambient |  |
| Salivary Cortisol Ambient |  |
| Thin Prep PAP Kit <i>Broom or Brush</i> Ambient |  |
| HPV Collection Device Ambient |  |
| ThinPrep Cytolyt Solution <i>Cyto FNA, rinse or brush</i> – Ambient |  |
| Formalin Containers 10 mL 10% Neutral Buffered Formalin 60 mL 10% Neutral Buffered Formalin 120 mL 10% Neutral Buffered Formalin Ambient |  |