

Validated Specimen Matrix

Collection Kit		Sample Type	Aptima Chlamydia trachomatis	Aptima Neisseria gonorrhoeae	Aptima Mycoplasma genitalium	Aptima Trichomonas vaginalis	Aptima Bacterial vaginosis	Aptima Candida vaginitis	Aptima HSV 1&2	Aptima hrHPV & Genotype
ThinPrep Pap Vial		Cervicovaginal	X	X		X				X
		Anal/rectal								X
Aptima Multitest Swab ¹		Vaginal (Clinician or Patient Collect ²)	X	X	X	X	X	X	X ³	
		Throat	X	X						
		Anal/rectal	X	X						
		Anogenital or Vesicular Lesion							X	
Aptima Unisex Swab		Endocervical	X	X	X	X				
		Male Urethral	X	X	X	X			X	
		Anogenital or Vesicular Lesion							X	
Aptima Urine Specimen		Female Urine	X	X	X	X			X	
		Male Urine	X	X	X	X			X	

¹If ordering testing for more than 3 infections, it is recommended that two swabs are collected and submitted.

²Vaginal specimens may only be self-collected within a healthcare setting where counseling is available to describe the self-collection procedure.

³Clinician-collect only and requires a separate swab

ACL Test Order Codes

Test Order Code	Directory Link	Test Name	Acceptable Collection			
			ThinPrep	Aptima MultiTest	Aptima Unisex	Aptima Urine
LAB10332 ¹ or Pap Order	https://apps.aurora.org/dos/public/lab/lab10332	ThinPrep Image-Guided Liquid Based Pap Test (GYN)	Cervicovaginal			
		ThinPrep Image-Guided Liquid Based Pap Test with reflex to Aptima HPV				
		ThinPrep Image-Guided Liquid Based Pap Test and Aptima HPV				
		ThinPrep Image-Guided Liquid Based Pap Test and Aptima HPV with reflex to HPV Genotype 16, 18/45				
LAB9939 or HPVHR	https://apps.aurora.org/dos/public/lab/lab9939	Aptima HPV	Cervicovaginal Anal/Rectal Vulvar			
LAB9938 or HPVGT	https://apps.aurora.org/dos/public/lab/lab9938	Aptima HPV Genotype 16, 18/45	Cervicovaginal Anal/Rectal Vulvar			
LAB9053 or CYTOLOGY	https://apps.aurora.org/dos/public/lab/cytology	Cytology Non-Gynecological	Anal/Rectal			

¹Chlamydia, Gonorrhea and Trichomonas can be added

Test Order Code	Directory Link	Test Name	Acceptable Collection			
			ThinPrep	Aptima MultiTest	Aptima Unisex	Aptima Urine
LAB9913 or CGRNA	https://apps.aurora.org/dos/public/lab/lab9913	Chlamydia trachomatis / Neisseria gonorrhoeae by NAA	Endocervical	Vaginal Throat Rectal	Endocervical Male Urethral	Urine
LAB9912 or CTRNA	https://apps.aurora.org/dos/public/lab/lab9912	Chlamydia trachomatis by NAA	Endocervical	Vaginal Throat Rectal	Endocervical Male Urethral	Urine
LAB9928 or NGRNA	https://apps.aurora.org/dos/public/lab/lab9928	Neisseria gonorrhoeae by NAA	Endocervical	Vaginal Throat Rectal	Endocervical Male Urethral	Urine
LAB9964 or TVRNA	https://apps.aurora.org/dos/public/lab/lab9964	Trichomonas Vaginalis by NAA	Endocervical	Vaginal	Endocervical Male Urethral	Urine
LAB9960 or SWOMG	https://apps.aurora.org/dos/public/lab/swomg	SwabOne Mycoplasma genitalium by NAA		Vaginal	Endocervical Male Urethral	Urine
LAB11872	https://apps.aurora.org/dos/public/lab/lab11872	Sexually Transmitted Infections Panel by NAA (Chlamydia trachomatis, Neisseria gonorrhoeae, Trichomonas vaginalis - 3 results)	Endocervical	Vaginal	Endocervical Male Urethral	Urine
LAB9957 or SWOBV	https://apps.aurora.org/dos/public/lab/swobv	SwabOne Vaginitis by NAA (Lactobacillus species, Gardnerella vaginalis and Atopobium vaginae - 1 test result)		Vaginal		
LAB9958 or SWOCN	https://apps.aurora.org/dos/public/lab/swocn	SwabOne Candida/Trichomonas Panel by NAA (Candida species, C. glabrata, Trichomoniasis - 3 test results)		Vaginal		
LAB9961 or SWOPNL	https://apps.aurora.org/dos/public/lab/swopnl	SwabOne Vaginitis Panel by NAA (Bacterial vaginosis, Candida vaginitis, Trichomoniasis - 4 test results)		Vaginal		
LAB9934 or HSVPCR	https://apps.aurora.org/dos/public/lab/lab9934	Herpes Simplex Virus 1 & 2 by NAA		Anogenital or Vesicular lesion Vaginal ²	Anogenital or Vesicular lesion Male Urethral	Urine

²Clinician-collect only