

Collection Kit		Sample Type	Aptima CT/NG	Aptima Mycoplasma genitalium	Aptima Trichomonas vaginalis	Aptima HSV 1&2	Aptima BV	Aptima CV/TV	Aptima HPV	Aptima HPV Genotype
ThinPrep Pap Vial		Cervicovaginal	X		X				X	X
		Anal/rectal							X	X
Aptima Multitest Swab		Vaginal (Clinician and Patient Collect)	X	X	X	X	X	X		
		Throat	X							
		Anal/rectal	X							
		Anogenital Lesion				X				
Aptima Unisex Swab		Endocervical	X	X	X					
		Male Urethral	X	X	X	X				
		Anogenital Lesion				X				
Aptima Urine Specimen		Female Urine	X	X	X	X				
		Male Urine	X	X	X	X				

Test Order Code	Directory Link	Test Name	Acceptable Collection				Sample Type
			ThinPrep	Aptima MultiTest	Aptima Unisex	Aptima Urine	
LAB10332* or Pap Order	<a href="https://apps.aurora.org/dos/public/lab/lab10332">https://apps.aurora.org/dos/public/lab/lab10332</a>	ThinPrep Image-Guided Liquid Based Pap Test (GYN)	X				Cervicovaginal
		ThinPrep Image-Guided Liquid Based Pap Test with reflex to Aptima HPV					
		ThinPrep Image-Guided Liquid Based Pap Test and Aptima HPV					
		ThinPrep Image-Guided Liquid Based Pap Test and Aptima HPV with reflex to HPV Genotype 16, 18/45					
LAB9939 or HPVHR	<a href="https://apps.aurora.org/dos/public/lab/lab9939">https://apps.aurora.org/dos/public/lab/lab9939</a>	Aptima HPV	X				Cervicovaginal, anal/rectal, vulvar
LAB9938 or HPVGT	<a href="https://apps.aurora.org/dos/public/lab/lab9938">https://apps.aurora.org/dos/public/lab/lab9938</a>	Aptima HPV Genotype 16, 18/45	X				Cervicovaginal, anal/rectal, vulvar
LAB9053 or CYTOLOGY	<a href="https://apps.aurora.org/dos/public/lab/cytology">https://apps.aurora.org/dos/public/lab/cytology</a>	Cytology Non-Gynecological	X				Anal/rectal, vulvar

\*Note Chlamydia, Gonorrhea and Trichomonas can be added

Test Order Code	Directory Link	Test Name	Acceptable Collection				Sample Type
			ThinPrep	Aptima MultiTest	Aptima Unisex	Aptima Urine	
LAB9913 or CGRNA	<a href="https://apps.aurora.org/dos/public/lab/lab9913">https://apps.aurora.org/dos/public/lab/lab9913</a>	Chlamydia trachomatis / Neisseria gonorrhoeae by NAA	X	X	X	X	Vaginal, Endocervical, Male Urethral, Throat, Rectal, Urine
LAB9912 or CTRNA	<a href="https://apps.aurora.org/dos/public/lab/lab9912">https://apps.aurora.org/dos/public/lab/lab9912</a>	Chlamydia trachomatis by NAA	X	X	X	X	Vaginal, Endocervical, Male Urethral, Throat, Rectal, Urine
LAB9928 or NGRNA	<a href="https://apps.aurora.org/dos/public/lab/lab9928">https://apps.aurora.org/dos/public/lab/lab9928</a>	Neisseria gonorrhoeae by NAA	X	X	X	X	Vaginal, Endocervical, Male Urethral, Throat, Rectal, Urine
LAB9964 or TVRNA	<a href="https://apps.aurora.org/dos/public/lab/lab9964">https://apps.aurora.org/dos/public/lab/lab9964</a>	Trichomonas Vaginalis by NAA	X	X	X	X	Vaginal, Endocervical, Male Urethral, Urine
LAB9960 or SWOMG	<a href="https://apps.aurora.org/dos/public/lab/swomg">https://apps.aurora.org/dos/public/lab/swomg</a>	SwabOne Mycoplasma genitalium by NAA		X	X	X	Vaginal, Endocervical, Male Urethral, Urine
LAB9957 or SWOBV	<a href="https://apps.aurora.org/dos/public/lab/swobv">https://apps.aurora.org/dos/public/lab/swobv</a>	SwabOne Bacterial Vaginosis by NAA (Lactobacillus species, Gardnerella vaginalis and Atopobium vaginae - 1 test result)		X			Vaginal
LAB9958 or SWOCN	<a href="https://apps.aurora.org/dos/public/lab/swocn">https://apps.aurora.org/dos/public/lab/swocn</a>	SwabOne Candida/Trichomonas Panel by NAA (Candida species, C. glabrata, Trichomoniasis - 3 test results)		X			Vaginal
LAB9961 or SWOPNL	<a href="https://apps.aurora.org/dos/public/lab/swopnl">https://apps.aurora.org/dos/public/lab/swopnl</a>	SwabOne Vaginitis Panel by NAA (Bacterial vaginosis, Candida vaginitis, Trichomoniasis, Mycoplasma genitalium - 5 test results)		X			Vaginal
LAB9934 or HSVPCR	<a href="https://apps.aurora.org/dos/public/lab/lab9934">https://apps.aurora.org/dos/public/lab/lab9934</a>	Herpes Simplex Virus 1 & 2 by NAA		X	X	X	Anogenital lesion, Vaginal, Male Urethral, Urine